

RollOn manual transfer aid

Support while standing



With a Safe Working Load of 160 kg, RollOn is a versatile sit-to-stand manual assistive device, supporting patients in mobilisation to and from a seated position, such as transfers to and from the toilet, or between a bed and chair.

Features and benefits

Padded rear supports provide security and are easily positioned and adjusted

Safety lock ensures rear supports remain in an ideal position

Width-adjustable base for optimal adaptability

Crafted from aluminium, which is both **lightweight and environmentally-friendly**

Lockable dual castors ensure stable positioning and usage

RollOn manual transfer aid



Adjustable lower leg support



Manual base-width adjustment



Safety lock



Adjustable handle/seat frame

Technical Specifications

Maximum Load	160 kg
Weight	25 kg
Height	1410 mm
Base	
height	110 mm
length	890 mm
width	620 - 900 mm
width expansion	manual
Turning circle	1050 mm
Castors	
front	75 mm double band
rear	75 mm double band with brakes
Material	aluminium
Product number	6001001

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

RollOn is a manual assistive device developed for safe transfers to and from the sitting position. Class I
 Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
 Doc. No: 5EN183319-01, 25 October 2016

Enhancing outcomes for patients and their caregivers:

Hill-Rom